

## **Isolation of DNA from Filters using Hot Detergent and CTAB**

2X Lysis buffer is 100mM Tris pH 8.0, 40mM EDTA, 100mM NaCl, 1% SDS. Preheat to 70C.

Add zircon beads, about 50 $\mu$ L, to the filter along with 100 $\mu$ L lysis buffer. (You can adjust the amount of lysis buffer depending upon your sample volume – you want to cover your whole sample with buffer. When working with tissue pieces, you can also use a small eppendorf mortar and pestle to grind the tissues, or make your own by melting the tips of 1000 $\mu$ L pipet tips.)

Vortex and then heat @ 70C for 5 minutes.

Vortex sample for 30 seconds on the homogenize setting.

Heat @ 70C for 5 minutes.

Vortex as before.

Transfer the solution to a new microfuge tube and make the solution at least 0.7M with NaCl (if 100 $\mu$ L, add 35 $\mu$ L 3M NaCl). Add 10% CTAB to make a 1% solution (15 $\mu$ L) and incubate at 60 – 70C for 10 minutes.

Add an equal volume of phenol/chloroform (200 $\mu$ L), vortex, then spin @ 12,000 rpm for 5 minutes.

Remove the aqueous layer to a new eppendorf tube.

Precipitate supernatants using 100% isopropanol (0.6 volume isopropanol, no need to add extra salt.) Chill for at least 1 hour @ -80C, or overnight at -20C, before spinning @ 12,000 rpm for 10 minutes. Check bottom of tube carefully for very small whitish pellet. Air dry for 15 minutes.

Resuspend the pellet in 10 – 20 $\mu$ L sterile milliQ water or 0.01M Tris pH 7.6, and determine concentration by A260 if desired.