









Switch Accessible iPad Apps
Dan Herlihy – dherlihy@juno.com

The following apps of the iTunes app store are switch accessible. Listed is the name, cost, which interface will work with it and a screen shot of the product.


<p>Scene & Heard (49.99) Switch accessible with R.J. Coopers Bluetooth Switch Interface</p> <p>Change switch access settings in the Settings panel of the iPad</p>	
<p>Predictable (159.99) Switch accessible with R.J. Coopers Bluetooth Switch Interface</p> <p>Change switch access settings in the Settings panel of the iPad</p>	
<p>TapSpeak Choice (149.99) Switch accessible with AbleNet Blue2 Bluetooth switch or R.J. Coopers Bluetooth Switch Interface</p> <p>Change switch access settings in the Settings panel of the iPad</p>	
<p>TapSpeak Button (29.99) Switch accessible with R.J. Coopers Bluetooth Switch Interface</p> <p>Change switch access settings in the application</p>	
<p>SoundingBoard (49.99) Switch accessible with AbleNet Blue2 Bluetooth switch</p> <p>Change switch access settings in the application</p>	

<p>Rad Sounds – (29.99) Switch access using R.J. Coopers Blue Tooth Switch Interface</p> <p>Change switch access settings in the application</p>	
<p>TapToTalk (subscription based) Switch access using R.J. Coopers Blue Tooth Switch Interface</p> <p>Change switch access settings in the Settings panel of the iPad</p>	
<p>Tap Speak Sequence – (29.99 – 49.99 for Plus version) Switch access using R.J. Coopers Bluetooth Switch Interface or AbleNet Blue2 Bluetooth switch</p> <p>Change switch access settings in the application</p>	
<p>GoTalk Now (79.99) Switch access using R.J. Coopers Blue Tooth Switch Interface or Attainment Switch App</p>	
<p>Alexicom AAC – to be used in conjunction with published Alexicom AAC pages http://www.alexicomaac.com/dnn/Default.aspx</p>	
<p>Five Sharks Swimming – Inclusive (2.99) Switch accessible with AbleNet Blue2 Bluetooth switch & R.J. Coopers Blue Tooth Switch Interface</p>	
<p>Inclusive Smarty Pants (2.99) Switch accessible with AbleNet Blue2 Bluetooth switch & R.J. Coopers Blue Tooth Switch Interface</p>	

<p>Aunty Maggie's Recipe – Inclusive (2.99) Switch accessible with AbleNet Blue2 Bluetooth switch & R.J. Coopers Blue Tooth Switch Interface</p>	
<p>Attainment Switch (4.99) Turns iPad into a wireless switch interface for Mac or PC computer</p>	

Where do you change the settings?





Before you begin, both switch interfaces (R.J. Coopers Bluetooth Switch Interface or AbleNet Blue2 Bluetooth switch) are Bluetooth devices Bluetooth needs to be turned on in the Settings panel of the iPad. Then, the switch interface must be turned on for the iPad to see it, and selected in the Bluetooth control panel for the switch to connect to the iPad.




Shown here is a screen shot from the Settings window of the iPad. Scene & Heard is the app that is selected from the list on the left, and the settings for it are shown in the window to the right. Options are available to choose the Access Method – Direct Touch or Switch Access, Scanning Rate, etc. Switch settings for Predictable, Tap Speak Choice and Tap to Talk are also found in the Settings panel.



Switch access settings for the other applications are listed here, and a screen shot of the menu to open for settings:

<p>Tap Speak Sequence</p>	
<p>Sounding Board</p>	
<p>Five Sharks Swimming, Auntie Maggie's Recipe and Smarty Pants all have a similar looking switch menu screen – the only difference is the name of the activity</p>	
<p>Tap Speak Button</p>	

Rad Sounds – switch



R.J Cooper Bluetooth Switch Interface (\$99.00)

- 2 switch ports send a ~1 command for port 1 and a ~3 command for port 2, built-in rechargeable battery,



Ablenet Blue2 Bluetooth Switch (\$149.00)

- 5 pre-sets on the right panel of the device also allow for Pg Up/Dn, Left/Right Arrow, Up/Down Arrow, Right/Left Mouse Clicks, Tab – Return
- 2 built in switches (pedals) send default Shift/Enter settings, 2 external switch ports
- Uses 2 AA batteries



it-iClick - iPad Switch Interface - <http://www.inclusive.co.uk/it-iclick-p5716>

it-iClick is a Simply Works interface specifically designed for connection to the iPad. It plugs directly into the Apple port and allows up to six Simply Works products to be paired with the iPad to achieve switch access and keyboard functions in single or multi-user configurations.

